






# How It Works

- A soft sleeve at the entry point of the pump forms an airtight seal at the base of penis.
- Pumping also puts pressure on the ligaments in the area, moderately increasing length.
- Pumping expels air from the cylinder to create a vacuum.
- This vacuum effect draws bloodflow into penis, causing it to swell temporarily.

**BEFORE USING YOUR PUMP :**  
Make sure to connect the pump hose to the pump cylinder. Locate the quick-release valve button just above the pressure gauge (Power Pump) or just above the squeeze trigger (Supreme Pump). If you feel you have applied too much pressure or experience discomfort at any time, press the quick-release valve to release pressure and withdraw your penis from the cylinder.



## Care and Maintenance

Disassemble, wash each piece by hand with soap and warm water, and air dry completely before reassembling and storing.

**WARNING:** If you have symptoms of erectile dysfunction (i.e., inability to achieve an erection that is sufficient for sexual intercourse), consult your physician prior to using this device to avoid a potentially harmful delay in diagnosing any of the most common causes of this condition, such as diabetes, multiple sclerosis, cirrhosis of the liver, chronic renal failure, or alcoholism. This pump is not intended for medical use. It has not been clinically tested nor is it approved. Use of this product is done at your own risk. This is a recreational device and is not empirically validated.

# The Benefits of Pumping

- Boost confidence
- Straighten curved shaft
- Increase stamina
- Larger, fuller penis head
- Moderately enhance length and girth
- Counteract erectile dysfunction
- Please partners for more satisfying sex and relationships
- Instant visible results that last long enough for sex
- Create rock-hard erections on demand

## Fast Facts

- The first pumps were medical devices for men living with erectile dysfunction, with the earliest models dating back to the 1870s.
- Pumping costs less and carries fewer risks and side effects in comparison to any other enlargement method.
- Using a c-ring when pumping can help you maintain a harder, stronger erection.

# Staying Safe

Every man is different, and you should always consult your doctor before beginning any new health and wellness regimen. However, following these standard recommendations will improve your chances of properly and safely attaining optimal results.

- Pay attention to the pressure gauge. Too much pressure for too long can have negative effects, so follow the instructions closely and stop or hit the quick-release button when you reach the recommended pressure limit. A safe range is between 0 in.Hg and 5 in.Hg.
- Don't overpump. It may seem logical to increase the pressure when you don't feel the pump working as intensely as you did at the start of your pumping session, but go by the pressure gauge reading. The body naturally adjusts to the effects of pumping, and it's normal for sensation to change the longer you pump.
- If you feel any pain or notice red spots (broken capillaries) appear, stop pumping and hit the quick-release button immediately.

OPTIMIZE  
EVERY  
OPPORTUNITY

DOCJOHNSON

OptiMALE

Men's Sexual Health  
and Wellness

Pumps

For more information, scan with  
your smartphone.



optimaleproducts.com

Copyright © 2017. Doc Johnson Enterprises.  
All Rights Reserved.

docjohnson.com

KIN-0692-95  
3-2017