

MAINSQUEEZE™

VARIABLE PRESSURE ULTRASKYN™ STROKER

DOCJOHNSON®

Girls of Social Media™



@playmateiryna
5210-01-BX



@brittanya187
5210-02-BX



@viking.barbie
5210-03-BX



@ana.lorde
5210-04-BX



@layna.me
5210-05-BX



Bailey Rayne
5230-01-BX

Building on the success of Doc Johnson's® Main Squeeze™ line of hard-case ULTRASKYN™ strokers featuring lifelike molds of top adult performers, the Main Squeeze™: Girls of Social Media collection brings consumers a new kind of star-studded cast straight from their favorite mainstream sites and apps.

Boasting followers in the multi-millions, the Girls of Social Media put an exciting new face on the design and functions that make Main Squeeze™ so popular: a built-in squeeze plate designed to vary the tightness of each stroke, unique internal textures in lifelike ULTRASKYN™, and a suction-adjusting end cap that twists to different levels of intensity.

Main Squeeze's™ sturdy case and screw-on top keep it clean, discreet and convenient for travel. All items in the Girls of Social Media collection include detailed instructions on usage and care.



Featured Products

@playmateiryna - 5210-01-BX

Molded directly from Playmate and Instagram star @playmateiryna.

@brittanya187 - 5210-02-BX

Molded directly from internationally-recognized Instagram star @brittanya187.

@viking.barbie - 5210-03-BX

Molded directly from cover model and Instagram star @viking.barbie.

@ana.lorde - 5210-04-BX

Molded directly from voluptuous, tattooed Instagram star @ana.lorde.

@layna.me - 5210-05-BX

Molded directly from model and instagram star @layna.me.

Bailey Rayne - 5230-01-BX

Molded directly from fan favorite cam star Bailey Rayne.

Care and Maintenance

Gently pull from the top to remove the ULTRASKYN™ insert from inside the case. Remove squeeze plate section from magnetic dock (or from inside the case if it falls in during disassembly.) Clean each section with Main Squeeze™ Toy Cleaner or mild soap and warm water. Allow all parts to air dry completely and dust the ULTRASKYN™ tunnel with Main Squeeze™ Refresh Powder.

Dimensions

	Total Length	Width/Diameter	Circumference
ALL ITEMS	8" (20.32 cm)	3" (7.62 cm)	9.5" (24.13 cm)

CONTACT INFO:

www.docjohnson.com